

Aflatoxins/Mycotoxins ELISA Kits

ELISA kits for the determination of afla/mycotoxins in food/feed/beverage samples



Some cereals, fruits, other vegetal products, milk, cheese and complex food, feed and beverages based on these products can hide residuals of aflatoxins and mycotoxins, whose ingestion can provoke cancer or other serious diseases or tissue damages to humans and animals.

The Aflatoxins/Mycotoxins ELISA Kits proposed by Astori Tecnica for the quantitative determination of these toxins are extremely easy to use. All of them are based on cutting edge technology, with excellent analytical performances such as high specificity and exquisite sensitivity; their reliability is great with a wide range of possible commodities, from grains and flours to complex food, feed and beverage samples.

Their 2 available sizes with individually breakable wells avoid wasting test material and money.

Advantages:

- ☺ Wide range of quantitative ELISA Kits, with excellent sensitivity and specificity
- ☺ Quick extraction times and swift test protocols
- ☺ Ready-to-use reagents
- ☺ Suitable to different kinds of grains, flours, food, feed, dairy and beverage samples
- ☺ Availability of kits for "difficult" matrices and for specific kinds of samples
- ☺ 48-well and 96-well formats, with individually breakable wells
- ☺ Long shelf-life

Product range:

Code No.	Product description	No. wells	Applications
941AFL01M96	Total Aflatoxins ELISA Kit	96	Grains, cereals, nuts, feed
941AFL01M48	Total Aflatoxins ELISA Kit	48	Grains, cereals, nuts, feed
981AFL01LM96	Total Aflatoxins (LOW MATRIX) ELISA Kit	96	Grains, cereals, silage, nuts, spices, feed, etc.
981AFL01LM48	Total Aflatoxins (LOW MATRIX) ELISA Kit	48	Grains, cereals, silage, nuts, spices, feed, etc.
941BAFL01B196	Aflatoxin B1 (Rapid) ELISA Kit	96	Grains, cereals, nuts, feed
941BAFL01B148	Aflatoxin B1 (Rapid) ELISA Kit	48	Grains, cereals, nuts, feed
981BAFL01LM96	Aflatoxin B1 (LOW MATRIX) ELISA Kit	96	Grains, cereals, silage, nuts, spices, feed, etc.
981BAFL01LM48	Aflatoxin B1 (LOW MATRIX) ELISA Kit	48	Grains, cereals, silage, nuts, spices, feed, etc.
961AFLM01M96	Aflatoxin M1 ELISA Kit	96	Milk and dairy products
991AFLM01U96	Aflatoxin M1 ELISA Kit, for urine	96	Urine
991AFLM01U48	Aflatoxin M1 ELISA Kit, for urine	48	Urine
941DON01M96	Deoxynivalenol (DON) ELISA Kit	96	Grains, cereals, feed
941DON01M48	Deoxynivalenol (DON) ELISA Kit	48	Grains, cereals, feed
941OCH01M96	Ochratoxin ELISA Kit	96	Grains, cereals, feed
941OCH01M48	Ochratoxin ELISA Kit	48	Grains, cereals, feed
981OCH01W48	Ochratoxin ELISA Kit, for wine, qualitative/semi-quantitative test	48	Wine, grapes, juice, grape must
981OCH01ALC96	Ochratoxin ELISA Kit, alcoholic beverages	96	Wine, grapes, juice, grape must, beer
981OCH01ALC48	Ochratoxin ELISA Kit, alcoholic beverages	48	Wine, grapes, juice, grape must, beer
961OCH01COF96	Ochratoxin A ELISA Kit, for coffee, cocoa and spices	96	Coffee, cocoa, cocoa butter, spices
961OCH01COF48	Ochratoxin A ELISA Kit, for coffee, cocoa and spices	48	Coffee, cocoa, cocoa butter, spices
991OCH01MS96	Ochratoxin A ELISA Kit, for milk and serum	96	Human and animal milk, serum and plasma
991OCH01MS48	Ochratoxin A ELISA Kit, for milk and serum	48	Human and animal milk, serum and plasma
951FUM01C96	Fumonisin ELISA Kit	96	Corn
951FUM01C48	Fumonisin ELISA Kit	48	Corn
951ZEA01N96	Zearalenone ELISA Kit	96	Cereal crops and bread
951ZEA01N48	Zearalenone ELISA Kit	48	Cereal crops and bread

Distributed by: